BL21-CodonPlus (DE3)-RIPL Chemically Competent Cell 产品说明书

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BL21-CodonPlus (DE3)-RIPL Competent Cell	100µl /支
pUC19 (control vector, 10pg/µl)	10µl
保存条件 (保质期):	-80°C (6个月)

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$F^- ompT hsdS(r_B^- m_B^-) dcm^+ Tet^R gal \lambda(DE3) endA Hte [argU proL Cam^R] [argU ileY leuW Strep/Spec^R]$

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● Sample Induction Protocol (for reference only)

1. Inoculate a single colony from a freshly streaked plate into 3ml of LB medium containing the appropriate antibiotic for the plasmid and host strain.
2. Incubate with shaking at 200 rpm at 37°C overnight.
3. Inoculate 50 ml of LB medium containing the appropriate antibiotic with 0.5 ml of the overnight culture prepared in step 2 (use the 500 ml triangular flask as the container would be better).
4. Incubate with shaking at 150 rpm at 37°C until the OD 600 reaches 0.5-0.8. (0.6 recommended; about 2.5h).
5. (Optional) Pipet 1ml of the cultures into clean microcentrifuge tubes and place the tubes on ice until needed for gel analysis or storage at -20°C. These will serve as the non-induced control samples.
6. Add IPTG to a final concentration of 1 mM. Optimal time for induction of the target protein may vary from 2-16 hours, depending on the protein.
7. Incubate with shaking at 120 rpm at 37°C for 2-4 hours. To determine the optimal time for induction of the target protein, it is recommended that a time course experiment be performed varying the induction from 2-16 hours.
8. Place the culture on ice for 10 minutes. Harvest cells by centrifugation at 5,000×g for 10 minutes at 4°C.
9. Remove the supernatant and store the cell pellet at -20°C (storage at lower temperatures is also acceptable).

IPTG 配制:

Prepare a 1 M solution of IPTG (Isopropyl-β-D-thiogalactoside; Isopropyl-β-D-thiogalactopyranoside) by dissolving 2.38 g of IPTG in dd water and adjust the final volume to 10 ml. Filter sterilize before use.

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