



400 Pinnacle Way, Ste 450  
Norcross, GA 30071  
USA  
Tel: 770-931-8400  
Fax: 770-931-0230

## RNase-Free DNase I Set

The E.Z.N.A. RNase-Free DNase I Set is optimized for use with E.Z.N.A.® Total RNA protocols. Normally DNase I digestion is not required for RNA purified with HiBind® RNA Mini Columns as our silica-based spin column technology efficiently removes the majority of DNA without enzymatic digestion. However, certain sensitive RNA applications may require further DNA removal.

**Activity:** 10,000 Kunitz units/mg

One Kunitz unit is defined as the amount of DNase I that causes an increase in  $A_{260}$  of 0.001 per minute per milliliter at 25 °C, pH 5.0, with highly polymerized DNA as the substrate (1).

**Reaction Time:** 15-20 minutes on column at 20-30 °C

**Concentration:** 20 Kunitz/ $\mu$ L

**Storage/Stability:** Store at -20°C. All components of the RNase-Free DNase I Set are stable for at least 24 months from the date of purchase when stored at -20°C.

**Format:** Lyophilized enzyme with RNase-free buffer and water

### Kit Contents

Product	Preps	Units of DNase I
E1091-00	5	180
E1091-01	50	1,500
E1091-02	200	6,000

## On-membrane DNase I Digestion Protocol

The following protocol is a short procedure for On-Membrane DNase I digestion. Please take a few minutes to read the user manual accompanying the E.Z.N.A.<sup>®</sup> RNA Kit thoroughly to become familiar with the protocol. Prepare all materials required before starting the RNA isolation procedure to minimize RNA degradation. Follow the standard E.Z.N.A.<sup>®</sup> RNA protocol until the optional step for on-membrane DNase I digestion.

1. For each HiBind<sup>®</sup> RNA Mini Column, prepare the DNase I stock solution as follows:

Buffer	Volume per Prep
E.Z.N.A. <sup>®</sup> DNase I Digestion Buffer	73.5 $\mu$ L
RNase-free DNase I (20 Kunitz/ $\mu$ L)	1.5 $\mu$ L
Total Volume	75 $\mu$ L

### Important Notes:

- DNase I is very sensitive and prone to physical denaturing. **Do not vortex the DNase I mixture.** Mix gently by inverting the tube.
  - Freshly prepare DNase I stock solution right before RNA isolation.
  - Standard DNase buffers are not compatible with on-membrane DNase I digestion. The use of other buffers may affect the binding of RNA to the HiBind<sup>®</sup> matrix and may reduce RNA yields and purity.
  - All steps must be carried out at room temperature. Work quickly, but carefully.
2. Insert the HiBind<sup>®</sup> RNA Mini Column containing the sample into a 2 mL Collection Tube.
  3. Add 1/2 volume RNA Wash Buffer I (compared to the standard E.Z.N.A.<sup>®</sup> RNA protocol) to the HiBind<sup>®</sup> RNA Mini Column.
  4. Centrifuge at 10,000 x g for 1 minute.
  5. Discard the filtrate and reuse the Collection Tube.
  6. Add 75  $\mu$ L DNase I digestion mixture directly onto the surface of the membrane of the HiBind<sup>®</sup> RNA Mini Column.  
**Note:** Pipet the DNase I directly onto the membrane. DNA digestion will not be complete if some of the mixture is retained on the wall of the HiBind<sup>®</sup> RNA Mini Column.
  7. Let sit at room temperature for 15 minutes.
  8. Add 1/2 volume RNA Wash Buffer I (compared to the standard E.Z.N.A.<sup>®</sup> RNA protocol) to the HiBind<sup>®</sup> RNA Mini Column.
  9. Continue to the RNA Wash Buffer II wash step in the standard E.Z.N.A.<sup>®</sup> RNA protocol.