碧云天生物技术/Beyotime Biotechnology  
订货热线：400-168-3301或800-8283301  
订货e-mail：order@beyotime.com  
技术咨询：info@beyotime.com  
网址：http://www.beyotime.com

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1. Wang M, Zhao XR, Wang P, Li L, Dai Y, Huang H, Lei P, Zhu HF, Shen GX. Glucose regulated proteins 78 protects insulinoma cells (NIT-1) from death induced by streptozotocin, cytokines or cytotoxic T lymphocytes. *Int J Biochem Cell Biol.* 2007;39(11):2076-82.
2. Zhong W, Zou G, Gu J, Zhang J. L-arginine attenuates high glucose-accelerated senescence in human umbilical vein endothelial cells. *Diabetes Res Clin Pract.* 2010;89(1):38-45.
3. Lin Q, Ding X, Qiu F, Song X, Fu G, Ji J. In situ endothelialization of intravascular stents coated with an anti-CD34 antibody functionalized heparin-collagen multilayer. *Biomaterials.* 2010;31(14):4017-25.
4. Ling Y, Ye X, Ji H, Zhang Y, Lai Y, Peng S, Tian J. Synthesis and evaluation of nitric oxide-releasing derivatives of farnesylthiosalicylic acid as anti-tumor agents. *Bioorg Med Chem.* 2010;18(10):3448-56.
5. Zhang XP, Xiong JY, Jiao YX, Wang G, Zuo ZY. Involvement of mitochondrial ATP-sensitive potassium channels in etomidate preconditioning-induced protection in human myeloid HL-60 cells. *Environmental Toxicology and Pharmacology.* 2010;29: 320-322.
6. Zou CG, Zhao YS, Gao SY, Li SD, Cao XZ, Zhang M, Zhang KQ. Homocysteine promotes proliferation and activation of microglia. *Neurobiol Aging.* 2010;31(12):2069-79.
7. Fu R, Chen Z, Wang Q, Guo Q, Xu J, Wu X. XJP-1, a novel ACEI, with anti-inflammatory properties in HUVECs. *Atherosclerosis.* 2011 Nov; 219(1):40-8.
8. Zhong X, Xiu LL, Wei GH, Liu YY, Su L, Cao XP, Li YB, Xiao HP. Bezafibrate enhances proliferation and differentiation of osteoblastic MC3T3-E1 cells via AMPK and eNOS activation. *Acta Pharmacol Sin.* 2011 May;32(5):591-600.
9. Liu Y, He J, Jiang L, Wu H, Xiao Y, Liu Y, Li G, Du Y, Liu C, Wan J. Nitric oxide production is associated with response to brown planthopper infestation in rice. *J Plant Physiol.* 2011 May 15;168(8):739-45.
10. Zhang H, Li W, Wang G, Su Y, Zhang C, Chen X, Xu Y, Li K. The distinct binding properties between avian/human influenza A virus NS1 and Postsynaptic density protein-95 (PSD-95), and inhibition of nitric oxide production. *Virol J.* 2011 Jun 13;8:298.
11. Jin Y, Zheng S, Fu Z. Embryonic exposure to cypermethrin induces apoptosis and immunotoxicity in zebrafish (Daniorerio). *Fish Shellfish Immunol.* 2011 Apr-May;30(4-5):1049-54.
12. Song H, Yin W, Zeng Q, Jia H, Lin L, Liu X, Mu L, Wang R. Hemokinins modulate endothelium function and promote angiogenesis through neurokinin-1receptor. *Int J Biochem Cell Biol.* 2012 Sep; 44(9):1410-21.
13. Hu SQ, Ye JS, Zong YY, Sun CC, Liu DH, Wu YP, Song T, Zhang GY. S-nitrosylation of mixed lineage kinase 3 contributes to its activation after cerebral ischemia. *J Biol Chem.* 2012 Jan 20;287(4):2364-77.
14. Peng YL, Liu YN, Liu L, Wang X, Jiang CL, Wang YX. Inducible nitric oxide synthase is involved in the modulation of depressive behaviors induced by unpredictable chronic mild stress. *J Neuroinflammation.* 2012 Jul 6;9:75..
15. Zhang N, Gong L, Zhang H, Cao C. High glucose-induced dysfunction of endothelial cells can be restored by HoxA9EC. *Ann Vasc Surg.* 2012 Oct;26(7):1002-10.
16. Liu L, Liu J, Wong WT, Tian XY, Lau CW, Wang YX, Xu G, Pu Y, Zhu Z, Xu A, Lam KS, Chen ZY, Ng CF, Yao X, Huang Y. Dipeptidyl peptidase 4 inhibitor sitagliptin protects endothelial function in hypertension through a glucagon-like peptide 1-dependent mechanism. *Hypertension.* 2012 Sep;60(3):833-41.
17. Jia Y, Zhao Z, Xu M, Zhao T, Qiu Y, Ooi Y, Yang B, Rong R, Zhu T. Prevention of renal ischemia-reperfusion injury by short hairpin RNA of endothelin A receptor in a rat model. *Exp Biol Med (Maywood).* 2012
- Aug 1;237(8):894-902.
18. Zhong XN, Wang HH, Lu ZQ, Dai YQ, Huang JH, Qiu W, Shu YQ, Xu W, Cheng C, Hu XQ. Effects of Naoxitong on atherosclerosis and inducible nitric oxide synthase expression in atherosclerotic rabbit. *Chin Med J (Engl).* 2013 Mar;126(6):1166-70.
19. Jiao K, Zhang J, Zhang M, Wei Y, Wu Y, Qiu ZY, He J, Cao Y, Hu J, Zhu H, Niu LN, Cao X, Yang K, Wang MQ. The identification of CD163 expressing phagocytic chondrocytes in joint cartilage and its novel scavenger role in cartilage degradation. *PLoS One.* 2013;8(1):e53312.
20. Jin LF, Wang YY, Zhang ZD, Yuan YM, Hu YR, Wei YF, Ni J. Cytotoxicity and genome-wide microarray analysis of intestinal smooth muscle cells in response to hexavalentchromium induction. *Dongwuxue Yanjiu.* 2013 Jun 8;34(3):E93-E100.
21. Jiang Q, Zhou Z, Wang L, Shi X, Wang J, Yue F, Yi Q, Yang C, Song L. The immunomodulation of inducible nitric oxide in scallop Chlamys farreri. *Fish Shellfish Immunol.* 2013 Jan;34(1):100-8.
22. Zhang J, Niu N, Li B, McNutt MA. Neuron-derived IgG protects neurons from complement-dependent cytotoxicity. *J Histochem Cytochem.* 2013 Dec;61(12):869-79.
23. Song Y, Zhao XP, Song K, Shang ZJ. Ephrin-A1 is up-regulated by hypoxia in cancer cells and promotes angiogenesis of HUVECs through a coordinated cross-talk with eNOS. *PLoS One.* 2013 Sep 9;8(9):e74464.
24. Jiang XS, Ni YQ, Liu TJ, Zhang M, Jiang R, Xu GZ. Generation and characterization of immortalized rat retinal microglial cell lines. *J Neurosci Res.* 2014 Apr;92(4):424-31.
25. Wang T, Sun P, Chen L, Huang Q, Chen K, Jia Q, Li Y, Wang H. Cinnamtannin D-1 protects pancreatic  $\beta$ -cells from palmitic acid-induced apoptosis by attenuating oxidative stress. *J Agric Food Chem.* 2014 Jun 4;62(22):5038-45.
26. Wang ZC, Liu H, Bai YG, Yu JW, Zhang HJ, Cheng YP, Bao JX, Ren XL, Ma HZ, Ma J. 28-Day hindlimb unweighting reduces expression of Rho kinase and inhibits its effects in femoral artery of rat. *J Physiol Biochem.* 2015 Jun;71(2):205-16.
27. Cao Q, Jiang Y, Shi J, Xu C, Liu X, Yang T, Fu P, Niu T. Artemisinin inhibits the proliferation, migration, and inflammatory reaction induced by tumor necrosis factor- $\alpha$  in vascular smooth muscle cells through nuclear factor kappa B pathway. *J Surg Res.* 2015 Apr;194(2):667-78.
28. Li YH, Xu Q, Xu WH, Guo XH, Zhang S, Chen YD. Mechanisms of protection against diabetes-induced impairment of endothelium-dependent vasorelaxation by Tanshinone IIA. *Biochim Biophys Acta.* 2015 Apr;1850(4):813-23.
29. Zhao B, Zhou B, Bao L, Yang Y, Guo K. Alpha-Tomatine Exhibits Anti-inflammatory Activity in Lipopolysaccharide-Activated Macrophages. *Inflammation.* 2015 Oct;38(5):1769-76.
30. Yuan C, Zhang H, Wang W, Li Y, Yan R, Xu L, Song X, Li X. Transmembrane protein 63A is a partner protein of Haemonchus contortus galectin in the regulation of goat peripheral blood mononuclear cells. *Parasit Vectors.* 2015 Apr 9;8:211.
31. Zhao Q, Zhang J, Wang H. PGC-1 $\alpha$  overexpression suppresses blood pressure elevation in DOCA-salt hypertensive mice. *Biosci Rep.* 2015 May 21;35(3).
32. Liu XB, Yang BX, Zhang L, Lu YZ, Gong MH, Tian JK. An in vivo and in vitro assessment of the anti-inflammatory, antinociceptive, and immunomodulatory activities of Clematis terniflora DC. extract, participation of aurantiamide acetate. *J Ethnopharmacol.* 2015 Jul 1; 169:287-94.
33. Sun X, Jiao X, Ma Y, Liu Y, Zhang L, He Y, Chen Y. Trimethylamine N-oxide induces inflammation and endothelial dysfunction in human umbilical veinendothelial cells via activating ROS-TXNIP-NLRP3

- inflammasome. *Biochem Biophys Res Commun.* 2016 Dec 2; 481(1-2):63-70.
34. Li P, Zhang L, Zhang M, Zhou C, Lin N. Uric acid enhances PKC-dependent eNOS phosphorylation and mediates cellular ER stress: A mechanismfor uric acid-induced endothelial dysfunction. *Int J Mol Med.* 2016 Apr;37(4):989-97.
  35. Linghu K, Lin D, Yang H, Xu Y, Zhang Y, Tao L, Chen Y, Shen X. Ameliorating effects of ,8-cineole on LPS-induced human umbilical vein endothelial cell injury bysuppressing NF- $\kappa$ B signaling in vitro. *Eur J Pharmacol.* 2016 Oct 15;789:195-201.
  36. Ren Y, Tao S, Zheng S, Zhao M, Zhu Y, Yang J, Wu Y. Salvianolic acid B improves vascular endothelial function in diabetic rats with blood glucose fluctuations viasuppression of endothelial cell apoptosis. *Eur J Pharmacol.* 2016 Nov 15;791:308-315.
  37. Gadahi JA, Yongqian B, Ehsan M, Zhang ZC, Wang S, Yan RF, Song XK, Xu LX, Li XR. Haemonchus contortus excretory and secretory proteins (HcESPs) suppress functions of goat PBMCs invitro. *Oncotarget.* 2016 Jun 14;7(24):35670-35679.
  38. Tian X, Zhao L, Song X, Yan Y, Liu N, Li T, Yan B, Liu B. HSP27 Inhibits Homocysteine-Induced Endothelial Apoptosis by Modulation of ROS Production andMitochondrial Caspase-Dependent Apoptotic Pathway. *Biomed Res Int.* 2016;2016:4847874.
  39. Yao J, Zou Z, Wang X, Ji X, Yang J. Pinoresinol Diglucoside Alleviates oxLDL-Induced Dysfunction in Human Umbilical Vein Endothelial Cells. *Evid Based Complement Alternat Med.* 2016;2016:3124519.
  40. Liu J, Zhang W, Shi H, Yang K, Wang G, Wang P, Ji J, Chu PK. In situ plasma fabrication of ceramic-like structure on polymeric implant with enhanced surface hardness,cytocompatibility and antibacterial capability. *J Biomed Mater Res A.* 2016 May;104(5):1102-12.
  41. Gadahi JA, Li B, Ehsan M, Wang S, Zhang Z, Wang Y, Hasan MW, Yan R, Song X, Xu L, Li X. Recombinant protein of Haemonchus contortus 24 kDa excretory/secretory protein (rHcES-24) modulate the immunefunctions of goat PBMCs in vitro. *Oncotarget.* 2016 Dec 20;7(51):83926-83937.
  42. Gadahi JA, Ehsan M, Wang S, Zhang Z, Wang Y, Yan R, Song X, Xu L, Li X. Recombinant protein of Haemonchus contortus 14-3-3 isoform 2 (rHcftt-2) decreased the production of IL-4and suppressed the proliferation of goat PBMCs in vitro. *Exp Parasitol.* 2016 Dec;171:57-66.
  43. Li J, Geng XY, Cong XL. Corrigendum to "PGC-1 $\alpha$  ameliorates AngiotensinII-induced eNOS dysfunction in human aortic endothelial cells" [Vasc. Pharmacol. 83 (2016) 90-97]. *Vascul Pharmacol.* 2016 Dec; 87:248.
  44. Li Y, Yuan C, Wang L, Lu M, Wang Y, Wen Y, Yan R, Xu L, Song X, Li X. Transmembrane protein 147 (TMEM147): another partner protein of Haemonchus contortus galectin on thegoat peripheral blood mononuclear cells (PBMC). *Parasit Vectors.* 2016 Jun 23;9(1):355.
  45. Gao Y, Wang W, Cao J, Wang F, Geng Y, Cao J, Xu X, Zhou J, Liu P, Zhang S. Upregulation of AUFI is involved in the proliferation of esophageal squamous cell carcinoma throughGCH1. *Int J Oncol.* 2016 Nov;49(5):2001-2010.
  46. Zhang ZB, Ruan CC, Chen DR, Zhang K, Yan C, Gao PJ. Activating transcription factor 3 SUMOylation is involved in angiotensin II-induced endothelial cellinflammation and dysfunction. *J Mol Cell Cardiol.* 2016 Mar;92:149-57.
  47. Zhou G, Shi QS, Huang XM, Xie XB. Comparison of transcriptomes of wild-type and isothiazolone-resistant *Pseudomonas aeruginosa* by usingRNA-seq. *Mol Biol Rep.* 2016 Jun;43(6):527-40.
  48. Lu J, Kan S, Zhao Y, Zhang W, Liu J. Novel naproxen/esomeprazole magnesium compound pellets based on acid-independent mechanism: invitro and in vivo evaluation. *Drug Dev Ind Pharm.* 2016 Sep;42(9):1495-503.
  49. Liu C, Hao Y, Huang J, Li H, Yang Z, Zeng Y, Liu J, Li R. Ghrelin accelerates wound healing in combined radiation and wound injury in mice. *Exp Dermatol.* 2017 Feb;26(2):186-193.
  50. Deng Z, Wang Z, Jin J, Wang Y, Bao N, Gao Q, Zhao J. SIRT1 protects osteoblasts against particle-induced inflammatory responses and apoptosis in asepticprosthesis loosening. *Acta Biomater.* 2017 Feb; 49:541-554.
  51. Guo Q, Xu L, Li H, Sun H, Liu J, Wu S, Zhou B. Progranulin causes adipose insulin resistance via increased autophagy resulting from activated oxidativestress and endoplasmic reticulum stress. *Lipids Health Dis.* 2017 Jan 31;16(1):25.
  52. Guo J, Wang Y, Wang J, Zhang J, Fang J. Radiosensitization of oral tongue squamous cell carcinoma by nanosecond pulsed electric fields (nsPEFs). *Bioelectrochemistry.* 2017 Feb;113:35-41.
  53. Wei SJ, Cheng L, Liang ES, Wang Q, Zhou SN, Xu H, Hui LH, Ge ZM, Zhang MX. Poly(ADP-ribose) polymerase 1 deficiency increases nitric oxide production and attenuates aorticatherogenesis through downregulation of arginase II. *Clin Exp Pharmacol Physiol.* 2017 Jan; 44(1):114-122.

Version 2016.12.08