碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-168-3301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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- Xiao-Long W, Chuan-Ping Y, Kai X, Ou-Jv Q. Selenoprotein W depletion in vitro might indicate that its main function is not as an antioxidative enzyme. *iochemistry (Mosc)*. 2010;75(2):201-7.
- Chen L, Sun BB, Wang T, Wang X, Li JQ, Wang HX, Zhang SF, Liu DS, Liu L, Xu D, Ou XM, Chen YJ, Yang T, Zhou H, Wen FQ. igarette smoke enhances {beta}-defensin 2 expression in rat airways via nuclear factor- $\{\kappa\}$ B activation. *Eur Respir J*. 2010 Sep;36(3):638-45.
- Zhang X, Zhao X, Ma Z. PYDDT, a novel phase 2 enzymes inducer, activates Keap1-Nrf2 pathway via depleting the cellularlevel of glutathione. *Toxicol Lett*. 2010 Nov 10;199(1):93-101.
- Wu X, Zhu Y, Yan H, Liu B, Li Y, Zhou Q, Xu K. Isothiocyanates induce oxidative stress and suppress the metastasis potential of human non-smallcell lung cancer cells. *BMC Cancer*. 2010 Jun 9;10:269.
- Wang J, Sun P, Bao Y, Liu J, An L. Cytotoxicity of single-walled carbon nanotubes on PC12 cells. *Toxicol In Vitro*. 2011 Feb;25(1):242-50.
- Zhang J, Zhao QZ, Duan GL, Huang YC. Influence of sulphur on arsenic accumulation and metabolism in rice seedlings. *Environmental and Experimental Botany*,2011Aug,72(1),34-40.
- Bao HG, Wang XY, Li JY, Wu CX. Comparison of effects of hypoxia on glutathione and activities of related enzymes in livers of Tibetchicken and Silky chicken. *Poult Sci*. 2011 Mar;90(3):648-52.
- Zhao G, Yao-Yue C, Qin GW, Guo LH. Luteolin from Purple Perilla mitigates ROS insult particularly in primary neurons. *Neurobiol Aging*. 2012 Jan;33(1):176-86.
- Chen T, Zhang L, Yue JQ, Lv ZQ, Xia W, Wan YJ, Li YY, Xu SQ. Prenatal PFOS exposure induces oxidative stress and apoptosis in the lung of rat off-spring. *Reprod Toxicol*. 2012 Jul;33(4):538-45.
- Xu HL, Yu XF, Qu SC, Qu XR, Jiang YF, Sui da Y. Juglone, from *Juglans mandshruica* Maxim, inhibits growth and induces apoptosis in humanleukemia cell HL-60 through a reactive oxygen species-dependent mechanism. *Food Chem Toxicol*. 2012 Mar;50(3-4):590-6.
- Fan S, Yu Y, Qi M, Sun Z, Li L, Yao G, Tashiro S, Onodera S, Ikejima T. P53-mediated GSH depletion enhanced the cytotoxicity of NO in silibinin-treated human cervicalcarcinoma HeLa cells. *Free Radic Res*. 2012 Sep;46(9):1082-92.
- Wang J, Sun P, Bao Y, Dou B, Song D, Li Y. Vitamin E renders protection to PC12 cells against oxidative damage and apoptosis induced bysingle-walled carbon nanotubes. *Toxicol In Vitro*. 2012 Feb; 26(1):32-41.
- Liu JH, Qi ZM, Huang Q, Wei XL, Ke ZG, Fang YS, Tian YC, Yu ZL. Study of energetic-particle-irradiation induced biological effect on*Rhizopus oryzae* through Synchrotron-FTIR micro-spectroscopy. *Journal of Molecular Structure*. 2012Jan,1031(16),1-8.

14. Lou H, Jing X, Ren D, Wei X, Zhang X. Eriodictyol protects against H<sub>2</sub>O<sub>2</sub>-induced neuron-like PC12 cell death through activation of Nrf2/ARE signaling pathway. *Neurochem Int.* 2012 Jul;61(2):251-7.
15. Ou XH, Li S, Wang ZB, Li M, Quan S, Xing F, Guo L, Chao SB, Chen Z, Liang XW, Hou Y, Schatten H, Sun QY. Maternal insulin resistance causes oxidative stress and mitochondrial dysfunction in mouse oocytes. *Hum Reprod.* 2012 Jul;27(7):2130-45.
16. Li DY, Xue MY, Geng ZR, Chen PY. The suppressive effects of Bursopentine (BP5) on oxidative stress and NF- $\kappa$ B activation in lipopolysaccharide-activated murine peritoneal macrophages. *Cell Physiol Biochem.* 2012;29(1-2):9-20.
17. Wang Y, Wang L, Zhu Z, Ma W, Lei C. The molecular characterization of antioxidant enzyme genes in *Helicoverpa armigera* adults and their involvement in response to ultraviolet-A stress. *J Insect Physiol.* 2012 Sep;58(9):1250-8.
18. Yang L, Qu M, Wang Y, Duan H, Chen P, Wang Y, Shi W, Danielson P, Zhou Q. Trichostatin A inhibits transforming growth factor- $\beta$ -induced reactive oxygen species accumulation and myofibroblast differentiation via enhanced NF-E2-related factor 2-antioxidant response element signaling. *Mol Pharmacol.* 2013 Mar;83(3):671-80.
19. Liu BN, Yan HQ, Wu X, Pan ZH, Zhu Y, Meng ZW, Zhou QH, Xu K. Apoptosis induced by benzyl isothiocyanate in gefitinib-resistant lung cancer cells is associated with Akt/MAPK pathways and generation of reactive oxygen species. *Cell Biochem Biophys.* 2013 May;66(1):81-92.
20. Lin M, Xia B, Yang M, Gao S, Huo Y, Lou G. Characterization and antitumor activities of a polysaccharide from the rhizoma of *Menispermum dauricum*. *Int J Biol Macromol.* 2013 Feb;53:72-6.
21. Zhou Q, Duan H, Wang Y, Qu M, Yang L, Xie L. ROCK inhibitor Y-27632 increases the cloning efficiency of limbal stem/progenitor cells by improving their adherence and ROS-scavenging capacity. *Tissue Eng Part C Methods.* 2013 Jul;19(7):531-7.
22. Ji X, Wang H, Zhu J, Zhu L, Pan H, Li W, Zhou Y, Cong Z, Yan F, Chen S. Knockdown of Nrf2 suppresses glioblastoma angiogenesis by inhibiting hypoxia-induced activation of HIF-1 $\alpha$ . *Int J Cancer.* 2014 Aug 1;135(3):574-84.
23. Liu X, Chen L, Feng B, Liu G. Reversing effect of sorcin in the drug resistance of human nasopharyngeal carcinoma. *Anat Rec (Hoboken).* 2014 Feb;297(2):215-21.
24. Gao Y, Li W, Liu X, Gao F, Zhao X. Reversing effect and mechanism of soluble resistance-related calcium-binding protein on multidrug resistance in human lung cancer A549/DDP cells. *Mol Med Rep.* 2015 Mar;11(3):2118-24.
25. Han S, Lv X, Wang Y, Gong H, Zhang C, Tong A, Zhang B, Yao H. Effect and mechanism of peroxisome proliferator-activated receptor- $\gamma$  on the drug resistance of the U-87 MG/CDDP human malignant glioma cell line. *Mol Med Rep.* 2015 Aug;12(2):2239-46.
26. Xie X, Zhao Y, Ma CY, Xu XM, Zhang YQ, Wang CG, Jin J, Shen X, Gao JL, Li N, Sun ZJ, Dong DL. Dimethyl fumarate induces necroptosis in colon cancer cells through GSH depletion/ROS increase/MAPKs activation pathway. *Br J Pharmacol.* 2015 Aug;172(15):3929-43.
27. Gao Y, Li W, Liu X, Gao F, Zhao X. Reversing effect and mechanism of soluble resistance-related calcium-binding protein on multidrug resistance in human lung cancer A549/DDP cells. *Mol Med Rep.* 2015 Mar;11(3):2118-24.
28. Yuan ZY, Hu YL, Gao JQ. Brain Localization and Neurotoxicity Evaluation of Polysorbate 80-Modified Chitosan Nanoparticles in Rats. *PLoS One.* 2015 Aug 6;10(8):e0134722.
29. Wang C, Wang G, Liu H, Hou YL. Protective effect of bioactive compounds from *Lonicera japonica* Thunb. against H<sub>2</sub>O<sub>2</sub>-induced cytotoxicity using neonatal rat cardiomyocytes. *Iran J Basic Med Sci.* 2016 Jan;19(1):97-105.
30. Liu X, Li H, Liu L, Lu Y, Gao Y, Geng P, Li X, Huang B, Zhang Y, Lu J. Ethylation of arginine by PRMT1 regulates Nrf2 transcriptional activity during the antioxidative response. *Biochim Biophys Acta.* 2016 Aug;1863(8):2093-103.
31. Zhang T, Zhang C, Fan X, Li R, Zhang J. Effect of C-type natriuretic peptide pretreatment on in vitro bovine oocyte maturation. *In Vitro Cell Dev Biol Anim.* 2016 Oct 19. [Epub ahead of print]
32. Jia Y, Wang HD, Wang Q, Ding H, Wu HM, Pan H. GSH depletion and consequent AKT inhibition contribute to the Nrf2 knockdown-induced decrease in proliferation in glioblastoma U251 cells. *Oncol Rep.* 2017 Apr;37(4):2252-2260.
33. Han X, Xue X, Zhao Y, Li Y, Liu W, Zhang J, Fan S. Rutin-Enriched Extract from *Coriandrum sativum* L. Ameliorates Ionizing Radiation-Induced Hematopoietic Injury. *Int J Mol Sci.* 2017 Apr 29;18(5): pii: E942.
34. Cai Y, Zhao B, Liang Q, Zhang Y, Cai J, Li G. The selective effect of glycyrrhizin and glycyrrhetic acid on topoisomerase II  $\alpha$  and apoptosis in combination with etoposide on triple negative breast cancer MDA-MB-231 cells. *Eur J Pharmacol.* 2017 Aug 15;809:87-97.
35. Cen J, Zheng B, Bai R, Zhang L, Zhang F, Zhang X. Triterpenoids from *Aglaia abbreviata* exert cytotoxicity and multidrug resistant reversal effect in MCF-7/ADM cells via reactive oxygen species induction and P-glycoprotein inhibition. *Oncotarget.* 2017 Apr 20;8(41):69465-69476.
36. Zhang R, Shao M, Han X, Wang C, Li Y, Hu B, Pang D, Xie Z. ATP synthesis in the energy metabolism pathway: a new perspective for manipulating CdSe quantum dots biosynthesized in *Saccharomyces cerevisiae*. *Int J Nanomedicine.* 2017 May 19;12:3865-3879.

Version 2019.04.16