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- Luo JG, Li L, Kong LY. Preparative separation of phenylpropanoid glycerides from the bulbs of *Lilium lancifolium* by high-speed counter-current chromatography and evaluation of their antioxidant activities. *Food Chemistry*. 2012 Apr;131(3):1056-62.
- Zhang D, Luo M, Wang W, Zhao C, Gu C, Zu Y, Fu Y, Yao X, Duan M. Variation of active constituents and antioxidant activity in *Pyrola* [*P. incarnata* Fisch.] from different sites in Northeast China *Food Chemistry*. 2013 Dec;141(3):2213-9.
- Zhao P, Duan L, Guo L, Dou LL, Dong X, Zhou P, Li P, Liu EH. Chemical and biological comparison of the fruit extracts of *Citrus wilsonii* Tanaka and *Citrus*

medica L. Food Chem. 2015 Apr 15;173:54-60.

4. Gu Y, Tang X, Xie L, Meng G, Ji Y. Aliskiren improves endothelium-dependent relaxation of thoracic aorta by activating PI3K/Akt/eNOS signal pathway in SHR. Clin Exp Pharmacol Physiol. 2016 Apr;43(4):450-8.
5. Chen Y, Meng G, Bai W, Ma Y, Xie L, Altaf N, Qian Y, Han Y, Ji Y. Aliskiren protects against myocardial ischaemia-reperfusion injury via an endothelial nitric oxide synthase-dependent manner. Clin Exp Pharmacol Physiol. 2017 Feb;44(2):266-274.
6. Li W, Li Y, Sun R, Zhou S, Li M, Feng M, Xie Y. Dual character of flavonoids in attenuating and aggravating ischemia-reperfusion-induced myocardial injury. Exp Ther Med. 2017 Aug;14(2):1307-1314.
7. Zhang HT, Zhang ZJ, Mo WC, Hu PD, Ding HM, Liu Y, Hua Q3, He RQ. Shielding of the geomagnetic field reduces hydrogen peroxide production in human neuroblastoma cell and inhibits the activity of CuZn superoxide dismutase. Protein Cell. 2017 Jul;8(7):527-537.
8. Huang C, Cao X, Chen X, Fu Y, Zhu Y, Chen Z, Luo Q, Li L, Song X, Jia R, Yin Z, Feng B, Zou Y. A pectic polysaccharide from Ligusticum chuanxiong promotes intestine antioxidant defense in aged mice. Carbohydr Polym. 2017 Oct 15;174:915-922.
9. Tian R, Zhang LN, Zhang TT, Pang HY, Chen LF, Shen ZJ, Liu Z, Fang Q, Zhang SY. Association Between Oxidative Stress and Peripheral Leukocyte Telomere Length in Patients with Premature Coronary Artery Disease. Med Sci Monitor. 2017 Sep 11;23:4382-4390.
10. Song Q, Liu L, Yu J, Zhang J, Xu M, Sun L, Luo H, Feng Z, Meng G. Dihydromyricetin attenuated Ang II induced cardiac fibroblasts proliferation related to inhibitory of oxidative stress. Eur J Pharmacol. 2017 Jul 15;807:159-167.
11. Xie X, Zheng Y, Liu X, Cheng C, Zhang X, Xia T, Yu S, Wang M. Antioxidant Activity of Chinese Shanxi Aged Vinegar and Its Correlation with Polyphenols and Flavonoids During the Brewing Process. J Food Sci. 2017 Oct;82(10):2479-2486.
12. Hou N, Cai B, Ou CW, Zhang ZH, Liu XW, Yuan M, Zhao GJ, Liu SM, Xiong LG, Luo JD, Luo CF, Chen MS. Puerarin-7-O-glucuronide, a water-soluble puerarin metabolite, prevents angiotensin II-induced cardiomyocyte hypertrophy by reducing oxidative stress. N-S Arch Pharmacol. 2017 May;390(5):535-545.
13. Wang F, Xiong ZY, Li P, Yang H, Gao W, Li HJ. From chemical consistency to effective consistency in precise quality discrimination of Sophora flower-bud and Sophora flower: Discovering efficacy-associated markers by fingerprint-activity relationship modeling. J Pharmaceut Biomed. 2017 Jan 5;132:7-16.
14. Ma TJ, Lan DH, He SZ, Ye Z, Li P, Zhai W, Chen WQ, Huang Y, Fu Y, Sun A, Wang YB, Ye Z, Li JL, Gao Y, Yan XL, Li ZH. Nrf2 protects human lens epithelial cells against H₂O₂-induced oxidative and ER stress: The ATF4 may be involved. Exp Eye Res. 2018 Feb 2;169:28-37.
15. Zhang S, Xu L, Liu YX, Fu HY, Xiao ZB, She YB. Zhang S, Xu L, Liu YX, Fu HY, Xiao ZB, She YB. Nat Prod Bioprospect. 2018 Feb 27.
16. Wan X, Ahmad H, Zhang L, Wang Z, Wang T. Dietary enzymatically treated Artemisia annua L. improves meat quality, antioxidant capacity and energy status of breast muscle in heat-stressed broilers. J Sci Food Agr. 2018 Jan 9.
17. He TB, Huang YP, Huang Y, Wang XJ, Hu JM, Sheng J. Structural elucidation and antioxidant activity of an arabinogalactan from the leaves of Moringa oleifera. Int J Biol Macromol. 2018 Jan 30;112:126-133.
18. Hua Q, Chen C, Tel Zur N, Wang H, Wu J, Chen J, Zhang Z, Zhao J, Hu G, Qin Y. Metabolomic characterization of pitaya fruit from three red-skinned cultivars with different pulp colors. Plant Physiol Bioch. 2018 Mar 2;126:117-125.
19. He TB, Huang YP, Huang Y, Wang XJ, Hu JM, Sheng J. Structural elucidation and antioxidant activity of an arabinogalactan from the leaves of Moringa oleifera. Int J Biol Macromol. 2018 Jun;112:126-133.
20. Zhao C, Xia T, Du P, Duan W, Zhang B, Zhang J, Zhu S, Zheng Y, Wang M, Yu Y. Chemical Composition and Antioxidant Characteristic of Traditional and Industrial Zhenjiang Aromatic Vinegars during the Aging Process. Molecules. 2018 Nov 12;23(11). pii: E2949.

Version 2020.05.19