碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-1683301 或 800-8283301  
订货 e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

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- Chen L, Gin KY, He Y. Effects of sulfate on microcystin production, photosynthesis, and oxidative stress in *Microcystis aeruginosa*. Environ Sci Pollut Res Int. 2016 Feb;23(4):3586-95.
- Gu Y, Tang X, Xie L, Meng G, Ji Y. Aliskiren improves endothelium-dependent relaxation of thoracic aorta by activating PI3K/Akt/eNOS signalpathway in SHR. Clin Exp Pharmacol Physiol. 2016 Apr;43(4):450-8.
- Chen Y, Meng G, Bai W, Ma Y, Xie L, Altaf N, Qian Y, Han Y, Ji Y. Aliskiren protects against myocardial ischaemia-reperfusion injury via an endothelial nitric oxide synthase-dependent manner. Clin Exp Pharmacol Physiol. 2017 Feb;44(2):266-274.
- Lian W, Jia H, Xu L, Zhou W, Kang, Liu A, Du G. Multi-Protection of DL0410 in Ameliorating Cognitive Defects in D-Galactose Induced Aging Mice. Front Aging Neurosci. 2017 Dec 8;9:409.
- Mu X, Liu X, Wang X, Dai H, Liu C. Growth of rutile TiO<sub>2</sub> nanorods in Ti and Cu ion sequentially implanted SiO<sub>2</sub> and the involved mechanisms. Nanotechnology. 2018 Jan 1;29(2):025601
- Hua Q, Chen C, Tel Zur N, Wang H, Wu J, Chen J, Zhang Z, Zhao J, Hu G, Qin Y. Metabolomic characterization of pitaya fruit from three red-skinned cultivars with different pulp colors. Plant Physiol Bioch. 2018 May;126:117-125.

7. Zhu X,Li D,Du Y,He W,Lu Y.DNA hypermethylation-mediated downregulation of antioxidant genes contributes to the early onset of cataracts in highly myopic eyes.*Redox Biol.* 2018 Oct;19:179-189.
8. Wang SQ,Yang XY,Cui SX,Gao ZH,Qu XJ.Heterozygous knockout insulin-like growth factor-1 receptor (IGF-1R) regulates mitochondrial functions and prevents colitis and colorectal cancer.*Free Radic Biol Med.* 2019 Apr;134:87-98.
9. Ren L,Han F,Xuan L,Lv Y,Gong L,Yan Y,Wan Z,Guo L,Liu H,Xu B,Sun Y,Yang S,Liu L.Clusterin ameliorates endothelial dysfunction in diabetes by suppressing mitochondrial fragmentation.*Free Radic Biol Med.* 2019 Dec;145:357-373.
10. Jiao Q,Zhang W,Jiang Y,Jiang L,Chen X,Liu B.Study on the Interactions Between Caffeoylquinic Acids With Bovine Serum Albumin: Spectroscopy, Antioxidant Activity, LC-MSn, and Molecular Docking Approach.*Front Chem.* 2019 Dec 6;7:840.
11. Long W,Mu X,Wang JY,Xu F,Yang J,Wang J,Sun S,Chen J,Sun YM,Wang H,Zhang XD.Dislocation Engineered PtPdMo Alloy With Enhanced Antioxidant Activity for Intestinal Injury.*Front Chem.* 2019 Nov 15;7:784.
12. Jiang J,Qi L,Lv Z,Jin S,Wei X,Shi FDietary Stevioside Supplementation Alleviates Lipopolysaccharide-Induced Intestinal Mucosal Damage through Anti-Inflammatory and Antioxidant Effects in Broiler Chickens.*Antioxidants (Basel).* 2019 Nov 21;8(12). pii: E575.

Version 2020.05.19