XL10-Gold Chemically Competent Cell 产品说明书

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XL10-Gold Competent Cell	100µl /支
pUC19 (control vector, 10pg/µl)	10µl
保存条件 (保质期):	-80°C (6个月)

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Tet^R Δ(*mcrA*)183 Δ(*mcrCB-hsdSMR-mrr*)173 *endA1 supE44 thi-1 recA1 gyrA96 relA1 lac Hte* [F' *proAB lacI*^qZΔM15 *Tn10* (Tet^R) *Amy Cam*^R]

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● Stratagene standard protocol

1. Pre-chill a 14-ml BD Falcon polypropylene round-bottom tubes on ice. Preheat NZY+ broth to 42°C.
2. Thaw the cells on ice. When thawed, gently mix and aliquot 100 µl of cells into the pre-chilled tubes.
3. Add 4 µl of the β-ME (β 巯基乙醇) to the aliquot of cells.
4. Swirl the tubes gently. Incubate the cells on ice for 10 minutes, swirling gently every 2 minutes.
5. Add 0.1-50 ng of the experimental DNA (or 2 µl of a ligation mixture) to the aliquot of cells.
6. Swirl the tubes gently, then incubate the tubes on ice for 30 minutes.
7. Heat-pulse the tubes in a 42°C water bath for 30 seconds. The duration of the heat pulse is critical.
8. Incubate the tubes on ice for 2 minutes.
9. Add 0.9 ml of preheated (42°C) NZY+ broth and incubate the tubes at 37°C for 1 hour with shaking at 225-250 rpm.
10. Plate ≤200 µl of the transformation mixture on LB agar plates containing the appropriate antibiotic (and containing IPTG and X-gal if color screening is desired).
11. Incubate the plates at 37°C overnight. If performing blue-white color screening, incubate the plates at 37°C for at least 17 hours to allow color development (color can be enhanced by subsequent incubation of the plates for 2 hours at 4°C).

● Preparation of Media and Reagents

NZY+ Broth (per Liter)

- 10 g of NZ amine (Casein hydrolysate)
- 5 g of Yeast extract
- 5 g of NaCl
- Add deionized H₂O to a final volume of 1 liter
- Adjust to pH 7.5 using NaOH and then autoclave
- Add the following filter-sterilized supplements prior to use:
 - 12.5 ml of 1 M MgCl₂
 - 12.5 ml of 1 M MgSO₄
 - 20 ml of 20% (w/v) glucose (or 10 ml of 2 M glucose)

LB Agar (per Liter)

- 10 g of NaCl
- 10 g of Tryptone
- 5 g of Yeast extract
- 20 g of Agar
- Add deionized H₂O to a final volume of 1 liter
- Adjust pH to 7.0 with 5 N NaOH and then autoclave
- Pour into petri dishes (~25 ml/100-mm plate)

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